

Breast Pump Comparison



Model	Finesse™	Calypso Essentials	Advanced Personal Double	Finesse with Dottie Tote™	Pump in Style Advanced® Starter Set	Pump in Style® Advanced	S2 Plus	S1 Plus
SKU	101A01	63.00.240	101035078	101A07	101035077	57063	MM011305	MM011091
Adjustable Controls	Yes, separate dials for multiphase pumping	Yes, digital	Yes, single speed/suction knob, let-down button	Yes, larger dials for multiphase pumping	Yes, single speed/suction knob, let-down button	Yes, single speed/suction knob, let-down button	Yes, digital, 2 phase cycling with let-down button	Yes, digital, 2 phase cycling with let-down button
BPA Free	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Total Bottles	2 bottles	2 bottles	2 bottles	6 bottles	2 bottles	4 bottles	2 bottles	2 bottles
Built-In Bottle Holders	No	Yes	No	Yes (bag)	No	Yes (bag)	No	No
Closed System*	Yes	Yes	No	Yes	No	No	Yes	Yes
Flange Sizes	25mm	26mm	24mm	25mm and 30.5mm	24mm and 27mm	24mm	24mm and 28mm	24mm and 28mm
Universal Kit	Yes	No	No	Yes	No	No	No	No
Battery Pack*	Yes (6 - AA batteries)	Yes (6 - AA batteries)	No	Yes (6 - AA batteries)	Yes (8 - AA batteries)	Yes (8 - AA batteries)	No	Yes, Rechargeable Battery
Warranty	2 years	400 hours or 1 year, whichever is greater	1 year limited on motor, 90 days on other parts	1 year	1 year limited on motor, 90 days on other parts	1 year limited on motor, 90 days on other parts	2 years	2 years
Weight	1 lb.	4 lbs.	4 lbs.	1 lb.	5 lbs.	9 lbs.	2.5 lbs.	3 lbs.
WHO Code Compliant	Yes	No	No	Yes	No	No	Yes	Yes
Extras	N/A	N/A	N/A	Dottie Tote, cooler, extra bottles, valves and flanges	Extra valves and flanges	On-the-go Tote with built-in pump, cooler, 2 extra bottles	Hospital grade strength. Includes nightlight and timer	Hospital grade strength. Includes nightlight and timer
Covered by insurance?	Yes	Yes	Yes	Exclusive companies No additional costs	Exclusive companies No additional costs	Exclusive companies Accessory costs apply	Exclusive companies No additional costs	Exclusive companies Accessory costs apply

*Batteries not included, sold separately

**Closed System provides a barrier to protect milk from mold, bacteria and viruses