

Breast Pump Comparison



Model	Finesse™	Calypso Essentials	Advanced Personal Double	Easy Express 2	Finesse with Dottie Tote™ ¹	Pump in Style Advanced® Starter Set ¹	9 Plus	S2 Plus ¹	Breeze	S1 Plus ¹	Pump in Style® Advanced ¹
SKU	101A01	63.00.240	101035078	PA200DA	101A07	101035077	9Plus	MM011305	PA209DA	MM011091	57063
Adjustable Controls	Separate dials for multiphase pumping	Digital	Single speed/suction knob, let-down button	4 speed modes, 7 suction levels	Larger dials for multiphase pumping	Single speed/suction knob, let-down button	Digital,	Digital, 2 phase cycling, let-down button	Digital, 3 phases, 10 levels	Digital, 2 phase cycling let-down button	Single speed/suction knob, let-down button
BPA Free	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Total Bottles	2 bottles	2 bottles	2 bottles	2 bottles	6 bottles	2 bottles	2 bottles	2 bottles	2 bottles	2 bottles	4 bottles
Built-In Bottle	No	Yes	No	No	Yes (bag)	No	No	No	Yes (bag)	No	Yes (bag)
Closed System*	Yes	Yes	No	Yes	Yes	No	Yes	Yes	Yes	Yes	No
Flange Sizes	25mm	26mm	24mm	24mm	25mm and 30.5mm	24mm and 27mm	24mm and 28mm	24mm and 28mm	24mm	24mm and 28mm	24mm
Universal Kit	Yes	No	No	No	Yes	No	No	No	No	No	No
Battery Pack*	(6) AA batteries	(6) AA batteries	No	(4) AA batteries	(6) AA batteries	(8) AA batteries	Rechargeable	No	Rechargeable	Rechargeable	(8) AA batteries
Warranty	2 years	400 hours or 1 year, whichever is	1 year limited on motor, 90 days on	1 year limited on motor	2 years	1 year limited on motor, 90 days on other parts	2 years	2 years	1 year limited on motor	2 years	1 year limited on motor, 90 days on
Weight	1 lb.	4 lbs.	4 lbs.	1 lb.	1 lb.	5 lbs.	½ lb	2.5 lbs.	11.6 ounces	3 lbs.	9 lbs.
WHO Code Compliant	Yes	No	No	Yes	Yes	No	Yes	Yes	Yes	Yes	No
Extras	N/A	N/A	N/A	Mesh drying bag, carrying tote	Dottie Tote, cooler, extra bottles, valves and flanges	Extra valves and flanges	Ultra-quiet, extra flanges, multi-volt power for use in any country	Hospital grade strength. Includes nightlight and timer	Silicone cushions, mesh drying bag, carrying tote	Hospital grade strength. Includes nightlight and timer	On-the-go Tote with built-in pump, cooler, 2 extra bottles
Covered by insurance?	Yes	Yes	Yes	Yes	Yes, additional costs may apply	Yes, additional costs may apply	Yes, additional costs may apply	Yes, additional costs may apply	Yes, additional costs apply	Yes, additional costs apply	Yes, additional costs may apply

*AA Batteries not included, sold separately. Rechargeable batteries are built-in, included at no additional cost.

**Closed System provides a barrier to protect milk from mold, bacteria and viruses